

Dialysis Centers – Europe's Hospitals

ST-ELISABETH SITE – avenue De Frélaan 206, 1180 Uccle (High care) – Tel. +32 2 614 21 71
ST-MICHEL SITE – rue de Linthoutstraat 150, 1040 Etterbeek (High & Low cares) – Tel. +32 2 614 35 70
BELLA VITA MEDICAL CENTER – allée André Delvaux 16, 1410 Waterloo (Low care) – Tél. +32 2 614 42 22

About Mr – Mrs

1. Personal details

1.1. GENERAL

BIRTHDATE

HOME ADDRESS + PHONE NUMBER

HOLIDAY ADDRESS + PHONE NUMBER

Phone

Phone

Mobile

Mobile

DRIVER/TRANSPORT DETAILS

PERSON OF CONTACT

NEFROLOGISTS :

- Dr Fr. Reginster
- Dr A. Broughton
- Dr G. Clerbaux
- Dr M.-L. Cambier
- Dr A. Dejardin
- Dr L. Smets
- Dr A. Van Audenhove

1.2. MEDICAL HISTORY

NEPHROPATHY

ASSOCIATED DISEASES

CURRENT MEDICAL ISSUES

ALLERGY

WAITING LIST FOR DECEASED DONOR TRANSPLANTATION

☐ YES ☐ No

MEDICATIONS DURING DIALYSIS SESSIONS

☐ Intravenous ESA

☐ EPREX ☐ RECORMON ☐ ARANESP ☐ MIRCERA

Dosis / Week

☐ IV Iron

☐ INJECTAFER 100mg/2ml ☐ Other:

Dosis

☐ Other medication

Dosis

PERSONAL MEDICATIONS

--

LABORATORY RESULTS

Blood group		MRSA carrier	
HBS antigen		HCV	
HBS antibody		HIV	
CPE		VRE	

2. Dialysis treatment

2.1. PATIENT DATA

PROGRAM			Hours
---------	----------------------	----------------------	-------

DRY WEIGHT		kg	INTER-DIALYTIC WEIGHT GAIN		kg
------------	----------------------	----	----------------------------	----------------------	----

BP BEFORE DIALYSIS		AFTER DIALYSIS	
--------------------	----------------------	----------------	----------------------

PROBLEM DURING SESSIONS

--

VASCULAR ACCESS

Type ☐ AV fistula

- | | | |
|-------------------|---------------------------------------|---|
| ▪ <i>side</i> | ☐ left | ☐ right |
| ▪ <i>place</i> | ☐ arm | ☐ forearm ☐ wrist |
| ▪ <i>material</i> | ☐ needle (15g) | ☐ catheter |
| ▪ <i>type</i> | ☐ sharp | ☐ blunt |
| ▪ <i>mode</i> | ☐ uniponction | ☐ biponction |

☐ Central venous catheter

- | | | |
|----------------|----------------------------------|--------------------------------------|
| ▪ <i>side</i> | ☐ left | ☐ right |
| ▪ <i>place</i> | ☐ jugular | ☐ sub-clavian |

Slot volume (pure heparin):	arterial side		ml
	venous side		ml

2.2. TECHNICAL DATA

DIALYSES TYPE

☐ HD ☐ HDF (hemodiafiltration) ☐ predilution ☐ postdilution

DIALYSATE COMPOSITION

Sodium (Na⁺) ☐ Mmol Potassium (K⁺) ☐ 0 Mmol
 ☐ 2 Mmol
Calcium (Ca⁺⁺) ☐ 1.25 Mmol ☐ 3 Mmol
 ☐ 1.50 Mmol ☐ 4 Mmol
 ☐ 1.75 Mmol
Glucose ☐ g/l Bicarbonate (HCO₃⁻) Mmol

ANTICOAGULATION

Type ☐ Low molecular weight heparin

Brand name

Dosis

☐ Unfractionated Heparin LEO 5000 UI / 1 ml

Loading dosis

UI

Maintenance dosis

UI/hour

Stop before end of session

☐ 0

☐ 20

☐ 30

☐ 40 minutes

FILTER

AVERAGE BLOOD FLOW

ml / minute

IN CASE OF

UF volume >

Na or UF profiles

UF volume >

Isolated UF (30 min / 1 lit)

Others